



Ergotherapeutic Seat Wellbeing for your back



» Preserves the natural pelvic and spinal position when sitting


» Significantly reduces the pressure on the intervertebral discs

» **TILT ANGLE** and height of the seat adjustable



Standard

- Standard saddle shape
- Particularly popular among people with hips of average width
- Ensures a healthy and comfortable posture
- For the elimination and prevention of back pain




Seat height ideally matched to the appropriate body size:

Height of gaslifts	S	M	L
Body size	up to 1.60 m	up to 1.90 m	over 1.90 m

(Please consider when placing an order.)

Cutaway

- Comfortable for people with narrow hips



With
NARROW
curvature!



Article(s)	Color	REF
Bestseller: Bambach Saddle Seat Standard	Black	
Bambach Cutaway, genuine leather	Black	B204051
Bambach Cutaway, med. grade vinyl	Black	B203001

Genuine leather



White #032	Ivory #031	Cream #030	Pink #039	Aztec #027	Sun #029	Caribbean #035	Grey #026	Mandarin #037	Apple #034
Heron #025	Poppy #036	Tan #028	Kingfisher #024	Damson #038	Atlantic Blue #023	Ink #022	Forest #033	Black #021	

Medical Grade Vinyl (soft, abrasion resistant, disinfectable)

White #012	Almond #011	Silver Pearl #008	Pink #014	Sunshine #009	Ruby #017	Lagoon #015	Stone #006	Grey #003	Baltic #002
Aubergine #013	Blue Heron #004	Black #001							



POSTURE matters:



THE PROBLEM


Conventional seats force the spine into an unhealthy "C" shape. This unnatural posture increases the pressure on your spinal discs to an alarming extent, which can result in many problems, mainly in combination with pain in the back, neck, arms, hands or shoulders.

THE SOLUTION

Thanks to its patented back curvature, the Bambach Saddle Seat is the only seat that ideally supports your pelvis, allowing your spine to keep its natural "S" shape. While this avoids a poor sitting posture, you benefit from a natural, healthy posture.

Small

- Particularly popular with people with a small frame




SMALL
saddle shape!



Large

- Particularly popular with people of a larger build



LARGE
saddle shape!



Article(s)	Color	REF
Bambach Small, genuine leather	Black	B205051
Bambach Large, genuine leather	Black	B207051

Please consider your body size when placing an order.
Height of gaslifts: S up to 1.60 m M up to 1.90 m L over 1.90 m



Bambach® Foot height adjuster


Seat particularly for operation theatres

- Easy seat height adjustment without using your hands
- Flexible and quick position change
- Ideal for surgery or prolonged treatment
- For large, medium and small gas lift heights

Article(s)	Content(s)	REF
Foot height adjuster large	1 pc.	BSTE-2015
Foot height adjuster medium	1 pc.	BSTE-2014
Foot height adjuster small	1 pc.	BSTE-2013

Additional accessories:

Leather Care Oil	50 ml	554 304
------------------	-------	---------





Bambach® Backrest Smart

- Fully disinfectable
- Height adjustable
- Dimensions 22 x 15 cm



Article(s)

Bambach Backrest Smart
Bambach Backrest Large

REF

BBKR-1001
BBKR-3051

Bambach® Armrest Smart

Relief during longer treatments

- Perfect for endodontics and dental surgery
- Relieves tension in the shoulder muscles
- Swivel-mounted armrest
- Length 55 cm
- Available in black vinyl only



Article(s)

Bambach Armrest Smart, mounted on the right
Bambach Armrest Smart, mounted on the left

REF

BARM-1901
BARM-3901

Bambach® Armrest

Relief during longer treatments

- Ideal for longer treatments such as endodontics and surgeries
- Relief for the shoulder muscles
- Can be swiveled and retrofitted
- In two versions:
- Armrest short (Length 30 cm)

Armrest long (Length 55 cm)



Bambach® Armrest Microscope

With three movable arm axes

- Armrest freely movable in width and length (free-swinging)
- Height adjustable (fixed)
- Relieves shoulder, neck and arms



Article(s)

Armrest long, Penuine leather black, 55 cm (right)
Armrest long, genuine leather black, 55 cm (left)
Armrest short, genuine leather black, 30 cm (right)
Armrest short, genuine leather black, 30 cm (left)

REF

BARM-1051
BARM-3051
BARM-2051
BARM-4051

Article(s)

Bambach Armrest Microscope

REF

BARM-6000