

### 303



	6
Handstück · Handpiece (HP)	
330 104 603391 ...	
<b>303.104. ...</b>	•

⊙<sub>max.</sub> 15000 min<sup>-1</sup>/rpm

Scheiben-, Polierer- und Bürstenträger, rostfreier Stahl  
Mandrel for discs, polishers and brushes, stainless steel

### 305



	6	6
Größe · Size	∅ 1/10 mm	
	050	080
Handstück · Handpiece (HP)		
330 104 604391 ...		
<b>305.104. ...</b>	○050	●080

● = ⊙<sub>max.</sub> 10000 min<sup>-1</sup>/rpm

○ = ⊙<sub>max.</sub> 15000 min<sup>-1</sup>/rpm

Träger für Scheiben, Polierer und Bürsten, rostfreier Stahl

Mandrel for discs, polishers and brushes, stainless steel

### 301 L



	6
Handstück · Handpiece (HP)	
330 104 610415 ...	
<b>301L.104. ...</b>	•

⊙<sub>max.</sub> 15000 min<sup>-1</sup>/rpm

### 329



	6
Handstück · Handpiece (HP)	
330 104 610417 ...	
<b>329.104. ...</b>	•

⊙<sub>max.</sub> 25000 min<sup>-1</sup>/rpm

Spindelträger für Polierer, rostfreier Stahl

Spindle-shaped mandrel for polishers, stainless steel

### 329 A



	6
Handstück · Handpiece (HP)	
<b>329A.104. ...</b>	•

⊙<sub>max.</sub> 15000 min<sup>-1</sup>/rpm



### 326



		6	6
Größe · Size	Ø 1/10 mm	020	030
Handstück · Handpiece (HP)			
330 104 609000 ...			
		<b>326.104. ...</b>	
		020	030

⊙<sub>max</sub> 15000 min<sup>-1</sup>/rpm  
Träger für Kauflächenpolierer, rostfreier Stahl  
Mandrel for occlusal polishers, stainless steel

### 310



		6
Handstück · Handpiece (HP)		
330 104 608000 ...		
		<b>310.104. ...</b>
		•

⊙<sub>max</sub> 30000 min<sup>-1</sup>/rpm  
Pop-on Träger für Polierscheiben/Finierscheiben,  
rostfreier Stahl  
Pop-on mandrel for polishing/finishing discs, stainless  
steel

636

### 327



		6
Handstück · Handpiece (HP)		
330 104 615421 ...		
		<b>327.104. ...</b>
		•

⊙<sub>max</sub> 15000 min<sup>-1</sup>/rpm  
Spezialträger, rostfreier Stahl  
Special mandrel, stainless steel

### 318



		6
Handstück · Handpiece (HP)		
330 104 623442 ...		
		<b>318.104. ...</b>
		•

⊙<sub>max</sub> 30000 min<sup>-1</sup>/rpm  
Papierstreifenträger, rostfreier Stahl  
Mandrel for sandpaper strips, stainless steel

### 314



		6
Handstück · Handpiece (HP)		
330 104 622444 ...		
		<b>314.104. ...</b>
		•

⊙<sub>max</sub> 10000 min<sup>-1</sup>/rpm  
Papierstreifenträger, rostfreier Stahl  
Mandrel for sandpaper strips, stainless steel

### 305 L



		6
Handstück · Handpiece (HP)		
330 104 604395 ...		
		<b>305L.104. ...</b>
		•

⊙<sub>max</sub> 15000 min<sup>-1</sup>/rpm  
Träger mit Linksgewinde, rostfreier Stahl  
Mandrel with left-hand thread, stainless steel

● 329 L



6

Handstück · Handpiece (HP)



330 104 6104 18 ...

● 329L.104. ...

$n_{max}$  25000 min<sup>-1</sup>/rpm

Träger mit Linksgewinde, rostfreier Stahl

Mandrel with left-hand thread, stainless steel