



# Etch-Rite™

38% phosphoric acid etch gel



## Product Overview

Etch-Rite is a soft, thixotropic blue gel with handling characteristics most preferred by clinicians. It provides the optimal etch pattern on dentin and enamel surfaces to ensure mechanical retention of bonding agents, restorative resins, and resin cements.

- Dispenses through small-gauge needles
- Stays where placed
- Washes off with ease



**THE INDUSTRY STANDARD FOR DECADES**

*More than 20,000,000 applications every year*

*Available in a wide variety of packaging options*

|                |                                       |               |   |
|----------------|---------------------------------------|---------------|---|
| <b>ETCH</b>    | 4 x 1.2 mL / 1.6 gm syringes + 8 tips | <b>ET-50</b>  | 2 x 25 mL syringes (65 gm), + 5 x 3 mL empty syringes + 50 tips |
| <b>ET-6</b>    | 6 ml / 8.3 gm syringe                 | <b>ET-50R</b> | 2 x 25 mL syringes (65 gm)                                      |
| <b>ET-12</b>   | 12 mL / 16.6 gm syringe               | <b>25B20</b>  | Light blue, 25 ga. x 1/2", pre-bent tips, pkg. 20               |
| <b>ET-24</b>   | 24 x 1.2 mL syringes                  | <b>25B50</b>  | Light blue, 25 ga. x 1/2", pre-bent tips, pkg. 50               |
| <b>ET-TWIN</b> | 2 x 3 mL syringes + 25 tips           | <b>25B100</b> | Light blue, 25 ga. x 1/2", pre-bent tips, pkg. 100              |

# Etch-Royale™

37% phosphoric acid etch gel



## Product Overview

For clinicians who prefer a creamier gel that readily settles into dentin and enamel but does not run, Etch-Royale is the perfect choice. The darker blue color is easier to see in thin applications. Etch-Royale has all the same features as Etch-Rite, but the consistency is slightly creamier than its famous sister product.

- Creamy, thixotropic consistency
- Dark blue color
- Washes off easily



|             |  |              |  |
|-------------|--|--------------|--|
| <b>ER</b>   | 4 x 1.2 mL / 1.6 gm syringes + 20 tips | <b>ER50</b>  | 2 x 25 mL syringes (65 gm) + 5 x 3 mL empty syringes + 50 tips |
| <b>ER24</b> | 24 x 1.2 mL syringes                   | <b>ER50R</b> | 2 x 25 mL syringes (65 gm)                                     |