

ACE™ German Steel

Premium Instruments

by Modern Orthodontics

High Quality Surgical Grade German Steel

Designed for maximum comfort and safety

The quality of steel has a crucial effect on the quality of instruments.



German Steel is rust proof and provides durability and long lasting performance

New Pleasant glassy mat sand finish

All instruments are guaranteed against manufacturing defects

Proper handling and the correct use of cleaning solutions are the most important things to consider to ensure a long lifetime of the instrument





Weingart pliers slim
Cat.No. AC-101



Jaraback Pliers
Cat.No. AC-144



Weingart pliers universal
Cat.No. AC-102




Crown and Band
Contouring Pliers
Cat.No. AC-143



Posterior Band Remover
Cat.No.AC-130



Distal end cutter long
with safety hold
Cat.No. AC-110



Lingual Weingart Plier fine
Long Handle 60° Angled
Cat.No. AC-103



Distal end cutter mini
with safety hold
Cat.No. AC-111



Omega Loop Forming
Tweed Style
Cat.No. AC-140





Distal end cutter
Safety Hold
Long Handle
Cat.No. AC-112



Distal end cutter
Flush cut, safety Hold
Long Handle
Cat.No. AC-113




How Pliers Curved
Cat.No. AC-142



Pin & Ligature Cutter
TC
Cat.No. AC-120




Pin & Ligature Cutter
TC Angled 15°
Cat.No. AC-123



Pin & Ligature Cutter
Mini TC
Cat.No. AC-121



Reversed Ligature Cutter
Angled 40° With T.C
Cat. No. AC-124



Hard Wire Cutter TC
Cat.No. AC-122




Tweed Arch Forming
(0.9 mm thick)
Cat.No. AC-135



Step Bending Plier
Cat.No. AC-138



Hook Crimping Plier
Cat.No. AC-136



Nance Loop Forming Plier
Cat.No. AC-139



Separator Placing Plier
Cat.No. AC-137



Bow Pliers Straight
Cat.No. AC-141



Heavy Wire Cutter with sleeve
for extra oral use Cat.No. AC-242



Light Wire Bird
Beak Plier
Cat.No. AC-150



Adam Plier
Cat.No. AC-162




Light Wire Bird
Beak Long
Cat.No. AC-151



Young Universal Plier
with grooves
Cat.No. AC-163




Bracket Remover Straight
Cat.No.AC-131



Hammer Head Plier
(For Distal Cinch Back)
Cat.No. AC-160



Bracket Remover Curved
Cat.No. AC-132



Three Prong Plier
Cat.No. AC-164



Stop "V" Bending Plier
Cat.No. AC-161
For "V" stop in SS wire




Hollow Chop
Contouring Plier
Cat.No. AC-165



De La Rosa Plier
With Grooves
Cat.No. AC-166

Heavy Duty Plier for TPA,
lingual arch and extraoral
appliances
Cat.No. AC-152





Ligature Forming Plier
Cat.No.AC-174



Mosquito Forceps Curved
Cat.No. AC-169



Torquing Pliers set of two
Cat.No.AC-172



Mosquito Forceps Straight
Cat.No. AC-167



Mathieu Needle Holder
Narrow Tip
Cat.No.AC-212



Band Pusher and scaler instrument
Cat.No.AC-200



Torquing Plier with Key
Cat.No.AC-171



Optical Pliers
Cat.No.AC-115



Double-Ended Distal End Bender
Cat.No.AC-203



Ligature director and tucker
Cat.No.AC-201



Dontrix Gauge for Tension and compression force measurement of upto 16 Oz Cat.No..... AC-254



Separating Strip Holder
Cat.No.AC-226

Ideal for cutting stainless steel



Turret 0 Deg Torque with 7 grooves from 0.016"-0.021" rectangular wires
Cat.No. AC-251

Band Cutting Scissors.
bands and ligatures.
Cat.No.AC-223



Boone Gauge Stainless Steel
Cat.No.AC-222



Turret 5 Deg & 7 deg Torque with 5 grooves. 7 deg for 0.016"-0.017" wires. 5 deg for 0.018", 0.019" & .021" rectangular wires
Cat.No. AC-252



Plier Stand SS Foldable
Cat.No.AC-243



Boone Gauge SS with Lead
Cat.No.AC-250



Impression Trays
Cat.No.AC-227



MBT Bracket Height Gauge Aluminum
for .018" slot.....Cat.No.AC-221
for .022" slot.....Cat.No.AC-220




MBT Bracket Height Gauge
Movable Head.....AC-225



Bracket Bonding Tweezer with Slot Positioner
Cat.No.AC-231



Bracket Bonding Tweezer Curved for Molar Tubes
Cat.No.AC-232



Bite Stick High Heat
Cat.No.AC-240



Moulded High Heat Plastic
With Stainless Tips

Safe For Dry Sterilization



Fine Serrated Tips for
Reduced Risk of Slippage



Mershon Band Pusher
Cat.No.AC-241

Lightweight Hollow handle



Boley Gauge Stainless Steel
Cat.No.AC-224



Bracket Positive Positioner
Cat.No.AC-230

Precise placement of direct bond tubes or attachments made simple with the spring-loaded mechanism and small, grooved tips. The Positive Positioner is made of lightweight, anodized aluminum.

Access first and second molars easily
Slides open smoothly to release tube

