

9.24 Auxiliary Products and Accessories

Auxiliary Products

AlastiK™ Force Products

AlastiK™ Easy-To-Tie Ligatures

AlastiK™ Force Products

- **Not Made With Natural Rubber Latex**

- Exciting colors

- All colors are fade resistant for a minimum of one month in artificial saliva, 8% ethanol in water, cottonseed oil and carbonated citrated sodas.

AlastiK™ Easy-To-Tie Ligatures

AlastiK Easy-To-Tie Ligatures are designed with a 45° bend for easy and efficient tie-wing hook up. The unique angled shape reduces the range of movement needed for bracket ligation, improving/decreasing ligation time.



AlastiK Easy-To-Tie Ligature

Choose from 20 fade-resistant colors. One size fits most twin brackets.

Quantity 100 sticks per pack (1,000 ligatures).

Clear	406-870
Gray	406-871
Smoke	406-872
Light Blue	406-873
Dusky Blue	406-874
Teal	406-876
Lavender	406-877
Navy Blue	406-878
Kelly Green	406-880
Yellow	406-881
Deep Black	406-882
Silver	406-884
Green Metallic	406-887
Radical Red	406-979
Screaming Pink	406-980
Plum Crazy	406-981
Orange Blast	406-982
Extreme Green	406-983
Mediterranean Blue	406-985
Obscure	406-986
Mixed Assorted	406-888




A Modules

Mini-StiK™ Patient Sticks

A Modules

• Not Made With Natural Rubber Latex

- Injection molded
- Easy to use
- Comfortable to wear
- Packing options:
 - AlastiK™ MiKro-StiK™ System for unit-of-use dispensing
 - Mini-StiK™ patient sticks
 - QuiK-StiK™ with central ring-grip
 - Loose

			
Outer Diameter	A-1	A-2	A-3
Inches	.125	.094	.110
(mm)	(3.2 mm)	(2.4 mm)	(2.8 mm)
(Bracket)	(Full)	(Single)	(Mini)



Mini-StiK Patient Sticks

• Not Made With Natural Rubber Latex

- A-1 or A-3 modules
- 24 modules per stik
- Injection molded for smooth, comfortable surface
- Exciting colors

417 StiKs per bag
(10,008) A-1

Gray	406-391
Clear	406-392



42 StiKs per bag (1,008)

Colors	A-1	A-3
Gray	406-381	406-351
Clear	406-382	406-352
Smoke	406-903	406-943
Light Blue	406-906	406-946
Dusky Blue	406-907	406-947
Teal	406-909	406-949
Lavender	406-910	406-950
Kelly Green	406-916	406-956
Yellow	406-918	406-958
Silver	406-922	406-962
Green Metallic	406-935	406-975

Colors	A-1	A-3
Navy Blue	406-937	406-977
Deep Black	406-938	406-978
Assorted	406-924	-
Radical Red	406-221	406-232
Screaming Pink	406-222	406-233
Plum Crazy	406-223	406-234
Orange Blast	406-224	406-235
Extreme Green	406-225	406-236
Obscure	406-226	406-237
Mediterranean Blue	406-227	406-238

9.26 Auxiliary Products and Accessories

QuiK-StiK™

Loose Modules

B Modules

QuiK-StiK

- **Not Made With Natural Rubber Latex**

- A-1 or A-2 modules
- 60 modules per stick
- Central thumb ring
- A-1 is .125 in. (3.2 mm) outer diameter
- A-2 is .094 in. (2.4 mm) outer diameter

167 StiKs per bag
(10,020) A-1

Gray	406-476
Clear	406-477



17 StiKs per bag (1,020)

Colors	A-1	A-2
Gray	406-407	406-408
Clear	406-417	-
Smoke	406-703	-
Light Blue	406-706	-
Dusky Blue	406-707	-
Teal	406-709	-
Lavender	406-710	-
Kelly Green	406-716	-
Yellow	406-718	-
Silver	406-722	-
Green Metallic	406-735	-

Colors	A-1	A-2
Navy Blue	406-737	-
Deep Black	406-738	-
Assorted	406-724	-
Radical Red	406-241	-
Screaming Pink	406-242	-
Plum Crazy	406-243	-
Orange Blast	406-244	-
Extreme Green	406-245	-
Obscure	406-246	-
Mediterranean Blue	406-247	-

Loose A-1 Modules – 1,000 per package

- **Not Made With Natural Rubber Latex**

Colors	Force	A-1
Clear	Medium	406-011
Gray	Medium	406-007
Clear	Heavy	406-009

B Modules – 100 per package

- **Not Made With Natural Rubber Latex**

- Injection molded
- 4 sizes



Colors	Force	Center to Center Inches (mm)	B-1 .145 (3.7)	B-2 .175 (4.4)	B-3 .230 (5.8)	B-4 .325 (8.3)
Clear	Medium		406-021	406-022	406-023	406-027
Gray	Medium		406-091	406-092	406-093	406-094
Clear	Heavy		406-024	-	-	-

See Archwire Products section for coated and uncoated ligature wires and auxiliary ligature tie hooks.

C Modules Precut Chain Reel

C Modules

- **Not Made With Natural Rubber Latex**
- Three configurations
- Refer to each product for center to center measurement


C-1 Short Chain



***C-2** Long Chain



CK Continuous Chain



Packaging options include: bobbins, spools, and pre-cut lengths.

*C-2 only available in Precut C Module Chains

Precut C Module Chains (Injection molded)

12 module unit-of-use sections (10/pk)

- **Not Made With Natural Rubber Latex**

Colors	Force	Center to Center Inches (mm)	C-1	C-2	CK
			Mandibular .150 (3.8)	Maxillary .170 (4.3)	.115 (2.9)
Clear	Medium		406-031	406-032	406-035
Gray	Medium		406-029	406-030	406-036

Colors	Force	Center to Center Inches (mm)	.150	.175
			(3.8)	(4.4)
Clear	Heavy		406-033	406-034

Continuous C Module Chain Reel

- **Not Made With Natural Rubber Latex**

- 3 ft. (91 cm) per reel
- Die-cut chain
- Medium force



Colors	Force	Center to Center Inches (mm)	C-1	CK	
			.140 (3.6)	.110 (2.8)	
Clear	Medium		406-137	406-869	Package of 4 reels
Gray	Medium		406-037	406-669	Package of 4 reels
Clear	Medium		406-140	-	Package of 40 reels
Gray	Medium		406-040	-	Package of 40 reels

9.28 Auxiliary Products and Accessories

Continuous Chains Bobbin

Elastomeric Ligature with Guard Continuous Chains Bobbin

- **Not Made With Natural Rubber Latex**
- 15 ft. (457 cm) per bobbin

- Die-cut chain
- Color coordinated bobbin for easy identification

- Medium force
- 1 bobbin/pk



Colors	in. (mm)	Distance Center to Center	
		C-1 .140 (3.6)	CK .110 (2.8)
Gray		406-611	406-621
Clear		406-612	406-622
Green Metallic		406-655	406-645
Hot Red		406-656	406-646
Navy Blue		406-657	406-647
Deep Black		406-658	406-648
Smoke		406-673	406-623
Light Blue		406-676	406-626
Dusky Blue		406-677	406-627
Teal		406-679	406-629
Lavender		406-680	406-630
Kelly Green		406-686	406-636

Colors	in. (mm)	Distance Center to Center	
		C-1 .140 (3.6)	CK .110 (2.8)
Yellow		406-688	406-638
Silver		406-692	406-642
Radical Red		406-297	406-306
Screaming Pink		406-298	406-307
Plum Crazy		406-299	406-308
Orange Blast		406-300	406-309
Extreme Green		406-301	406-310
Obscure		406-302	406-312
Mediterranean Blue		406-303	406-313



AlastiK™ Bobbin Dispenser


- **Not Made With Natural Rubber Latex**
- Holds 16 bobbins
- Removable handle
- Replaceable blade

406-927

Replacement Blade Set

- Replacement Blade Set (2)
- 406-928

Elastomeric Ligature with Guard



• **Not Made With Natural Rubber Latex**

- .124 in. (3.1 mm) outer diameter
- Reduces tooth contact with brackets on mandibular teeth
- Injection molded memory material
- Disc dispenser, 10 modules per disc, 10 discs per package
- 406-430 Clear
- 406-429 Silver

Separator Modules

S Modules (Separators)



• **Not Made With Natural Rubber Latex**

- Pre-banding separators
- Three sizes
- Packaging options:
 - Loose
 - QuiK-StiK™

- Thermal sensitive
- Optimal effect when not over-stretched beyond:

S1	12 mm
S2	15 mm
S3	20 mm
SX	15 mm

Loose S Modules

• **Not Made With Natural Rubber Latex**

- Three sizes
- 100/pk*

Colors	Outer Diameter Inches (mm)	S1	S2	S3
		Anterior .140 (3.6)	Posterior .175 (4.4)	Posterior .210 (5.3)
Clear		406-081	406-082	406-085
Blue		-	406-084	-
Green		406-083	-	-

*Box available upon request




SX
Posterior
.190 in.
(4.8 mm)
Radiopaque
Dusky Blue
406-086 100/pk
406-087 1,000/pk

QuiK-StiK™ S Modules

• **Not Made With Natural Rubber Latex**

- 80 modules per stick, 13 sticks/pk
- Two sizes



Colors	Inches (mm)	Outer Diameter	
		S1 .145 (3.7)	S2 .175 (4.4)
Gray		406-411	406-410





9.30 Auxiliary Products and Accessories

K Modules

K Modules

- **Not Made With Natural Rubber Latex**

- For intermaxillary Class II and intramaxillary Class III applications
- 100/pk




Colors	Force	Center to Center Inches (mm)	 K-0	 K-1	 K-2	 K-3
			.400 (10.2)	.455 (11.6)	.650 (16.5)	.850 (21.6)
Clear	Medium	-	-	406-051	-	-
Gray	Medium	406-068	406-074	406-075	-	-
Clear	Heavy	-	-	406-054	-	406-056

KX Modules



- **Not Made With Natural Rubber Latex**

- Greater force delivery than K modules
- Designed for intermaxillary application
- Typical application is from a hook on distal molar to an archwire hook in the cuspid region
- 100/pk


Colors	Force	Center to Center Inches (mm)	 KX-1	 KX-2	 KX-3
			.235 (6.0)	.310 (7.9)	.415 (10.5)
Clear	Medium		406-060	406-057	406-062
Gray	Medium		406-064	406-059	406-065

Rotation Wedge

Action Line

EZ3 Holder

Rotation Wedge



Labial

• Not Made With Natural Rubber Latex

- Serves as wedge when archwire is ligated to opposing tie-wings
- Preformed tabs snap over tie-wings

Labial

100/pk

Gray

406-043

Lingual .010 in. (.25 mm) ligature wire

10/pk

Gray

390-150

Action Line Thread

• Not Made With Natural Rubber Latex

- For rotating and extruding teeth, closing minor space and other intraoral applications
- Clear polyurethane elastomeric thread
- Self-sealing bobbin

25 ft. (7.62 m) Bobbin

Diameter

.020 in. (.51 mm) **407-020**

.025 in. (.64 mm) **407-025**

.030 in. (.76 mm) **407-030**

EZ3 Holder



• Not Made With Natural Rubber Latex

- Pre-bonding lingual wire support
- Eliminates need for wax, ligature wire, floss threader or manual assistance
- 406-044 10/pk

Action Line Tubing

• Not Made With Natural Rubber Latex

- Hollow tubing for tighter knots
- Natural color elastomeric tubing
- Self-sealing bobbin

25 ft. (7.62 m) Bobbin

Diameter

.020 in. (.51 mm) **407-120**

.025 in. (.64 mm) **407-125**

.030 in. (.76 mm) **407-130**







9.32 Auxiliary Products and Accessories

People Pack Intraoral Elastics

People Pack Intraoral Elastics

- **Made With Natural Rubber Latex**
- **Coated With Cornstarch**
- Precision-cut, surgical-grade latex
- No color additives
- Size and force ratings are clearly marked on packaging

- Cross indexed by color and symbol
- Patient pack illustration allows written directions on surface

Sizes (I.D.) (Slack Diameters)			1/8 in. (3.2 mm)	3/16 in. (4.6 mm)	1/4 in. (6.4 mm)	5/16 in. (7.9 mm)	3/8 in. (9.5 mm)	5/8 in. (16 mm)	3/4 in. (19 mm)
			Pack Color						
Force Rating*		Pkg. Code	Purple	Green	Red	Blue	Orange	Aqua	Yellow
Light	2 oz. (56.7 g)		Carlos 404-116	Pete 404-126	Gary 404-136	Greg 404-146	Tom 404-156	Andy 404-166	Ernie 404-176
Medium-light	3 1/2 oz. (99.2 g)		Dave 404-216	Dwight 404-226	Chuck 404-236	Cliff 404-246	Will 404-256	Paul 404-266	Hal 404-276
Medium	4 oz. (113.4 g)		Cathy 404-516	Bill 404-526	Elliot 404-536	Louie 404-546	Joe 404-556	-	-
Medium-heavy	5 1/2 oz. (155.9 g)		-	Ken 404-426	-	-	-	-	-
Heavy	6 oz. (170.1 g)		Debbie 404-616	Bummer 404-626	Fred 404-636	Roger 404-646	George 404-656	-	-
Extra-heavy	8 oz. (226.8 g)		Gloria 404-716	Roberto 404-726	Fran 404-736	Anna 404-746	John 404-756	-	-

Packaging: All Intraoral sizes – 100 elastics per patient pack, 30 packs per box (Total 3,000 elastics)

*Force rating based upon stretching the elastic three times slack diameter