
At-Ease™


PASSIVE NITI CLIP

SELF-LIGATING


by Modern Orthodontics



- CNC Milled Slot
- NiTi Clip
- Metal Injection Molded
- 80 Mesh Pad



Modern Orthodontics



At-Ease™
PASSIVE NiTi CLIP
SELF-LIGATING SYSTEM

At-Ease™

Passive Self Ligating with NiTi Clip
Self Ligating Two- Piece Bracket
True Twin with four tie wings and
time tested 80 mesh

Metal Injection Moulded (MIM)


for intricate,
and smooth design.

Low Profile Design

for comfort and minimum occusal
interferences.

CNC Milled Slot

for precision
and smoothness




At-Ease™

PASSIVE NiTi CLIP
SELF-LIGATING SYSTEM

Advanced 17-4 low nickel alloy
for maximum strength in thinnest walls and zero slot collapse.

Exact in - out dimensions.

ZeroF ultra low friction slot
Advanced polishing process creates mirror finish by removing high surface peaks at molecular level.



Rounded tie wings
for patient comfort.

CNC milled slot
for precision and smoothness.

Centre scribe line
for accurate bracket placement in line with the long axis of tooth.




80 Gauge mesh pad
Compound contoured diagonal mesh pad for maximum bonding strength and zero rocking.

Palladium alloy braze
assures strongest possible brazed joint and resists electrolytic reaction that can cause brazing failure during use.

SecureLok Hooks
Offer best locking of chains and accessories. Due to thinnest hook stem design large amount undercut is present beneath the half round ball of hook and chains will not slip out in between the patient appointments.


Clip made of Nitinol
For effortless repeated opening & closing
Tested for more than 1000 cycles




At-Ease™

PASSIVE NiTi CLIP
SELF-LIGATING SYSTEM

We are sure that you won't miss your conventional brackets.....




Fully versatile | Use with or without 'O'ring | No training required




No special instruments to open or close | Use with or without 'O'ring




End to end rotation control



No Friction



Generous under tie wing area



At-Ease™

PASSIVE NiTi CLIP
SELF-LIGATING SYSTEM

MBT Prescription
Maxillary ▲

Tooth		Torque	Angle	Width	R/L	ID	.022
1|1	Central	+17°	+4°	.113"	R	●	BTA-M2-101
					L	●	BTA-M2-102
2|2	Lateral	+10°	+8°	.106"	R	●	BTA-M2-103
					L	●	BTA-M2-104
3|3	Cuspid with hook	0°	+9°	.119"	R	●	BTA-M2-105
					L	●	BTA-M2-106
3|3	Cuspid with hook	-7°	+9°	.119"	R	●	BTA-M2-109
					L	●	BTA-M2-110
4/5|4/5	Bicuspid with hook	-7°	0°	.110"	R	●	BTA-M2-113
					L	●	BTA-M2-114
4/5|4/5	Bicuspid	-7°	0°	.110"	R	●	BTA-M2-115R
					L	●	BTA-M2-115L

MBT Prescription
Mandibular ▼

Tooth		Torque	Angle	Width	R/L	ID	.022
½|½	Lower Anterior	-6°	0°	.093"	R	●	BTA-M2-119R
					L	●	BTA-M2-119L
3|3	Cuspid with hook	-6°	+3°	.108"	R	●	BTA-M2-120
					L	●	BTA-M2-121
3|3	Cuspid with hook	0°	+3°	.108"	R	●	BTA-M2-122
					L	●	BTA-M2-123
4|4	1st Bicuspid with hook	-12°	+2°	.110"	R	●	BTA-M2-128
					L	●	BTA-M2-129
4|4	1st Bicuspid	-12°	+2°	.110"	R	●	BTA-M2-130
					L	●	BTA-M2-131
5|5	2nd Bicuspid with hook	-17°	+2°	.110"	R	●	BTA-M2-136
					L	●	BTA-M2-137
5|5	2nd Bicuspid	-17°	+2°	.110"	R	●	BTA-M2-138
					L	●	BTA-M2-139

MBT Full Kits ▲▼

Kits	.022
Upper/Lower 5x5 Hook on 3*	BKA-M2-102*
Upper/Lower 5x5 Hooks on 3,4,5*	BKA-M2-104*

* Cuspid brackets in these Kits are having 0 deg torque in uppers and -6 deg torque in lowers.
For any other cuspid torque combination, please place order for individual brackets.