
d'art™

MINI CERAMIC

POLY CRYSTALLINE

by Modern Orthodontics




Low Profile Translucent

Dove Tail Grooves


Micro Crack Free

99.9% Pure Alumina





Modern Orthodontics



d'art™
MINI BRACKET
CERAMIC BRACKET SYSTEM

d'art™
Mini Size
CERAMIC BRACKET SYSTEM

Mini ceramic bracket system featuring small size and great strength

Mini Size


Smaller than most of ceramic brackets in the industry

Low Profile Design

for comfort and minimum occlusal interferences

Dove Tail Grooves

for pure mechanical bonding




d'art™
MINI BRACKET
CERAMIC BRACKET SYSTEM

Highest Quality Material
Made from high quality poly crystalline ceramic

Archwire Slot Design
Diamond and heat polished for superior sliding mechanics

Compound Contour Base
Provides simplified and accurate bracket placement

Dovetail Bonding Base
Provides superior mechanical bond strength and predictable debonding




Smooth, Strong and Contoured Tie Wings
For added patient comfort

Removable Color Coded ID Dot
For easy bracket selection

Ample Tie Wing Area
For easy bracket ligation

Ultra-Smooth, Strong Hooks
For added patient comfort and secure engagement of elastics







d'art™
MINI BRACKET
CERAMIC BRACKET SYSTEM







Modern Orthodontics




d'art™
MINI BRACKET
CERAMIC BRACKET SYSTEM

Bonds & De-Bonds Just Like A Metal Bracket



De-bonding tool to be used



Grasp the brackets on the
mesio-distal side close to the tooth




Wrong method of De-bonding



Twist the bracket 90° clockwise
to safely remove the bracket

Modern Orthodontics



d'art™ MINI BRACKET
CERAMIC BRACKET SYSTEM

Roth Prescription
Maxillary


	Tooth	Torque	Angle	Width	R/L	ID	.018	.022
1|1	Central	+12°	+5°	.127"	R	●	BTD-R8-001	BTD-R2-001
					L	●	BTD-R8-002	BTD-R2-002
2|2	Lateral	+8°	+9°	.125"	R	●	BTD-R8-003	BTD-R2-003
					L	●	BTD-R8-004	BTD-R2-004
3|3	Cuspid with Hook	0°	+11°	.128"	R	●	BTD-R8-005	BTD-R2-005
					L	●	BTD-R8-006	BTD-R2-006
4/5|4/5	1st & 2nd Bicuspid	-7°	0°	.128"	R	●	BTD-R8-009	BTD-R2-009
					L	●	BTD-R8-010	BTD-R2-010
4/5|4/5	1st & 2nd Bicuspid Hk	-7°	0°	.128"	R	●	BTD-R8-011	BTD-R2-011
					L	●	BTD-R8-012	BTD-R2-012

Roth Prescription
Mandibular

	Tooth	Torque	Angle	Width	R/L	ID	.018	.022
1/2|1/2	Lower Anterior	0°	0°	.110"	R	●	BTD-R8-014	BTD-R2-014
					L	●	BTD-R8-015	BTD-R2-015
3 | 3	Cuspid With hook	-11°	+5°	.126"	R	●	BTD-R8-016	BTD-R2-016
					L	●	BTD-R8-017	BTD-R2-017
4 | 4	1st Bicuspid	-17°	0°	.128"	R	●	BTD-R8-020	BTD-R2-020
					L	●	BTD-R8-021	BTD-R2-021
4 | 4	1st Bicuspid with hook	-17°	0°	.128"	R	●	BTD-R8-022	BTD-R2-022
					L	●	BTD-R8-023	BTD-R2-023
5 | 5	2nd Bicuspid	-22°	0°	.128"	R	●	BTD-R8-024	BTD-R2-024
					L	●	BTD-R8-025	BTD-R2-025
5 | 5	2nd Bicuspid with hook	-22°	0°	.128"	R	●	BTD-R8-026	BTD-R2-026
					L	●	BTD-R8-027	BTD-R2-027

Full Kits

Kits	.018	.022
Upper/Lower 5x5 Hook on 3	BKD-R8-102	BKD-R2-102
Upper/Lower 5x5 Hooks on 3,4 & 5	BKD-R8-104	BKD-R2-104



d'art™
MINI BRACKET
CERAMIC BRACKET SYSTEM

MBT Prescription
Maxillary ▲

	Tooth	Torque	Angle	Width	R/L	ID	.018	.022
1|1	Central	+17°	+4°	.125"	R	●	BTD-M8-101	BTD-M2-101
					L	●	BTD-M8-102	BTD-M2-102
2|2	Lateral	+10°	+8°	.122"	R	●	BTD-M8-103	BTD-M2-103
					L	●	BTD-M8-104	BTD-M2-104
3|3	Cuspid with hook	0°	+9°	.122"	R	●	BTD-M8-105	BTD-M2-105
					L	●	BTD-M8-106	BTD-M2-106
3|3	Cuspid with hook	-7°	+9°	.122"	R	●	BTD-M8-107	BTD-M2-107
					L	●	BTD-M8-108	BTD-M2-108
4/5|4/5	1st & 2nd Bicuspid	-7°	0°	.122"	R	●	BTD-M8-109	BTD-M2-109
					L	●	BTD-M8-110	BTD-M2-110
4/5|4/5	1st & 2nd Bicuspid Hk	-7°	0°	.122"	R	●	BTD-M8-111	BTD-M2-111
					L	●	BTD-M8-112	BTD-M2-112

MBT Prescription
Mandibular ▼

	Tooth	Torque	Angle	Width	R/L	ID	.018	.022
½|½	Lower Anterior	-6°	0°	.116"	R	●	BTD-M8-114	BTD-M2-114
					L	●	BTD-M8-115	BTD-M2-115
3|3	Cuspid with hook	0°	+3°	.122"	R	●	BTD-M8-116	BTD-M2-116
					L	●	BTD-M8-117	BTD-M2-117
3|3	Cuspid with hook	-6°	+3°	.122"	R	●	BTD-M8-118	BTD-M2-118
					L	●	BTD-M8-119	BTD-M2-119
4|4	1st Bicuspid	-12°	+2°	.122"	R	●	BTD-M8-120	BTD-M2-120
					L	●	BTD-M8-121	BTD-M2-121
4|4	1st Bicuspid with hook	-12°	+2°	.122"	R	●	BTD-M8-122	BTD-M2-122
					L	●	BTD-M8-123	BTD-M2-123
5|5	2nd Bicuspid	-17°	+2°	.122"	R	●	BTD-M8-124	BTD-M2-124
					L	●	BTD-M8-125	BTD-M2-125
5|5	2nd Bicuspid with hook	-17°	+2°	.122"	R	●	BTD-M8-126	BTD-M2-126
					L	●	BTD-M8-127	BTD-M2-127

Full Kits ▲▼

Kits	.018	.022
Upper/Lower 5x5 Hook on 3	BKD-M8-102*	BKD-M2-102*
Upper/Lower 5x5 Hooks on 3,4 & 5	BKD-M8-104*	BKD-M2-104*
Upper/Lower 5x5 Hook on 3	BKD-M8-202	BKD-M2-202
Upper/Lower 5x5 Hooks on 3,4 & 5	BKD-M8-204	BKD-M2-204

* Cuspid brackets in these Kits are having -7 deg torque in uppers and -6 deg torque in lowers.
Other kits have O deg torque in all cuspids.
For any other cuspid torque combination, please place order for individual brackets.

© Modern Orthodontics Doc#MO/CTG/F01 Issue #1 Rev.00 Dated 25/07/16